

ALL TEST™ SARS-CoV-2 och Influenza A+B Antigen Combo Rapid Test (Nasal Swab) Package Insert For self-testing

REF ISIN-525H English

A rapid test for the qualitative detection of SARS-CoV-2 Nucleocapsid Protein, Influenza A and Influenza B nucleoprotein antigens present in nasal swab specimen.

For *in vitro* diagnostic use.

【PROCEDURE】

Wash your hands with soap and water for at least 20 seconds before and after test. If soap and water are not available, use hand sanitizer with at least 60% alcohol.

Remove the cover of the tube with extraction buffer and place the tube in the tube holder in the box.

Nasal swab specimen Collection

1. Remove the sterile swab from the pouch. Do not touch the soft tip of the swab.

2. Insert the swab into your nostril until you feel slight resistance (Approx. 2cm up your nose). Slowly twist the swab, rubbing it along the insides of your nostril for 10-15 times against the nasal wall.

Note: This may feel uncomfortable. Do not insert the swab any deeper if you feel strong resistance or pain.

When the nasal mucosa is damaged or bleeding, nasal swab collection is not recommended.

If you are swabbing others, please wear a face mask. With children, you may not need to insert the swab so far into the nostril. For very young children you may need another person to steady the child's head while swabbing.

Gently remove the swab.

4. Using the same swab, repeat step 2 in your other nostril.

5. Withdraw the swab.

Specimen Preparation

1. Place the swab into the extraction tube, ensure it is touching the bottom and the swab to mix well. Press the swab head against to the tube and rotate the swab for 15-seconds.


2. Remove the swab while squeezing the swab head against the outside of the extraction tube.

3. Place the swab in a plastic bag.

4. Close the cap or fit the tube tip onto the tube.

Testing

1. Remove the test cassette from the sealed foil pouch and use it within one hour. Best results will be obtained if the test is performed immediately after opening the foil pouch.
2. Invert the specimen extraction tube and add 3 drops of extracted specimen to each sample well(S) of the test cassette and start timer. Do not move the test cassette during test development.
3. Read the result after 10 minutes. Do not read the result after 20 minutes.



Note: After test is completed, place all the components into plastic bag and tightly seal, then dispose according to local regulation.

【READING THE RESULTS】

Please carefully follow your local COVID guidelines/requirements.

POSITIVE SARS-CoV-2: Two colored lines appear in the COVID-19 window. One colored line should be in the control region (C) and another colored line should be in the test region (T).

POSITIVE Influenza A: Two colored lines appear in the FLU A+B window. One colored line should be in the control region (C) and another colored line should be in the Influenza A region (A).

POSITIVE Influenza B: Two colored lines appear in the FLU A+B window. One colored line should be in the control region (C) and another colored line should be in the Influenza B region (B).

POSITIVE Influenza A + Influenza B: Three colored lines appear in the FLU A+B window. One colored line should be in the control region (C) and two colored lines should be in the Influenza A region (A) and Influenza B region (B).

NOTE: The intensity of the color in the test line region (T/B/A) will vary based on the amount of SARS-CoV-2 and/or Influenza A+B antigen present in the sample. So dark a color of the test region (T/B/A) should be considered normal.

INTERFERING RESULTS: means it is very likely you have COVID-19 and/or Influenza A+B.

NEGATIVE: One colored line appears in the control region (C). No apparent colored line appears in the test line region (T/B/A).

Materials Provided

• Test cassette • Package insert • Sterile swab • Extraction tube (NaCl 5g/L, Tris 3g/L, Proclin 300 0.02%, BSA 5g/L, Triton-X100 2g/L, pH 8.5)

Materials Required But Not Provided

• Timer

【INTENDED USE】

The SARS-CoV-2 and Influenza A+B Antigen Combo Rapid

Testing

1. Remove the test cassette from the sealed foil pouch and use it within one hour. Best results will be obtained if the test is performed immediately after opening the foil pouch.
2. Place the test cassette on a flat and level surface.
3. Invert the specimen extraction tube and add 3 drops of extracted specimen to each sample well(S) of the test cassette and start timer. Do not move the test cassette during test development.
4. Read the result after 10 minutes. Do not read the result after 20 minutes.

You are unlikely to have COVID-19 and/or Influenza A/Influenza B. However, it is possible for this test to give a negative result that is incorrect (a false negative) in some people with COVID-19 and/or Influenza A/Influenza B. This means you could possibly still have COVID-19 and/or Influenza A/Influenza B even though the test is negative.

In addition, you can repeat the test after 1-2 days. If in case of suspicion, repeat the test after 1-2 days, as the coronavirus/influenza virus cannot be precisely detected in all phases of an infection.

Even with a negative test result, distance and hygiene must be observed, migration/traveling, attending events and etc. should follow your local COVID/influenza guidelines/requirements.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test.

【PRECAUTIONS】

Please review the information in this package insert before performing the test.

For self-testing *in vitro* diagnostic use only. Do not use after expiration date.

• Do not eat, drink or smoke in the area where the specimens or kits are handled.

• Do not touch the buffer in the kit. Carefully handle the buffer and do not touch skin or eyes, lines with plenty of liquid will cause immediate if contacting.

• Store in a dry place at 2-30 °C (36-86 °F), avoiding areas of excess moisture. If the foil packaging is damaged or has been opened, please do not use.

• This test kit is intended to be used as a preliminary test for self-testing.

Test for children should be under the guidance of an adult.

• Wash hands thoroughly before and after handling.

• Please ensure that an appropriate amount of samples are used for testing. Too much or too little sample size may lead to deviation of results.

【STORAGE AND STABILITY】

Store as packaged in the sealed pouch at room temperature. Do not refrigerate. Do not freeze. Do not store beyond the expiration date printed on the sealed pouch.

DO NOT FREEZE: Do not use beyond the expiration date.

【MATERIALS】

Materials Provided

• Test cassette • Package insert • Sterile swab

• Extraction tube (NaCl 5g/L, Tris 3g/L, Proclin 300 0.02%, BSA 5g/L, Triton-X100 2g/L, pH 8.5)

Materials Required But Not Provided

• Timer

【INTENDED USE】

The SARS-CoV-2 and Influenza A+B Antigen Combo Rapid

Test (Nasal Swab) is a single-use test kit intended to detect the SARS-CoV-2, Influenza A and Influenza B virus that causes COVID-19 and/or Influenza with self-collected nasal swab specimen. The test is intended for use in symptomatic respiratory specimens during the acute phase of infection. Positive results indicate the presence of SARS-CoV-2 and/or Influenza A/Influenza B antigens, but clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Positive results are indicative of the presence of SARS-CoV-2 and/or Influenza A+B. Individuals who test positive should self-isolate. Positive results do not rule out bacterial infection or co-infection with other viruses. Negative results do not preclude SARS-CoV-2 and/or Influenza A+B infection.

【SUMMARY】

The novel coronaviruses belong to the *beta* genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

Influenza (avian influenza) known as AI is a highly contagious acute infection in the respiratory tract. It is a communicable disease easily transmitted through the coughing and sneezing of aerosolized droplets containing live virus. Influenza outbreaks occur each year during the fall and winter months. Type A viruses are typically more prevalent than type B viruses and are associated with most serious seasonal epidemics, while type B infections are usually milder.

【PRINCIPLE】

The SARS-CoV-2 and Influenza A+B Antigen Combo Rapid Test (Nasal Swab) is a qualitative membrane-based immunoassay for the detection of SARS-CoV-2 nucleocapsid protein, Influenza A and Influenza B nucleoproteins antigens.

Positive results indicate the presence of SARS-CoV-2 and/or Influenza A/Influenza B antigens.

Negative results indicate the absence of SARS-CoV-2 and/or Influenza A/Influenza B antigens.

Test results are determined by the presence of specific antibodies to SARS-CoV-2 and/or Influenza A/Influenza B.

【LIMITATIONS】

Performance was evaluated with nasal swab specimens only, using the procedures provided in this package insert.

The SARS-CoV-2 and Influenza A+B Antigen Combo Rapid Test (Nasal Swab) will only indicate the presence of SARS-CoV-2 and/or Influenza A/Influenza B antigens in the specimen.

If the test result is negative, it is because the very early infection virus may not be detected. It is recommended to test again with a new kit. Negative results do not rule out SARS-CoV-2 and/or Influenza A/Influenza B infection in those who have been in contact with the virus.

Positive results of COVID-19 may be due to infection with non-SARS-CoV-2 coronavirus strains or other interference factors. A positive result for influenza A or B does not preclude an underlying co-infection with another pathogen, therefore the possibility of an underlying bacterial infection should be considered.

Failure to follow these procedures may alter test performance.

False negative results may occur if a specimen is

improperly collected or handled. False negative results may occur if inadequate levels of viruses are present in the specimen.

【PERFORMANCE CHARACTERISTICS】

The SARS-CoV-2 and Influenza A+B Antigen Combo Rapid Test (Nasal Swab) has been evaluated with specimens obtained from the patient. RT-PCR is used as the reference method for SARS-CoV-2 and Influenza A+B Antigen Combo Rapid Test (Nasal Swab).

Results are considered positive if RT-PCR indicated a positive result. Specimens were considered negative if RT-PCR indicated a negative result.

【SARS-CoV-2 Test】

Description

Test Level

SARS-CoV-2 and Influenza A+B Antigen Combo Rapid Test	RT-PCR (nasopharyngeal swab)	Total
Positive	161	2
Negative	5	482
Total	166	484

96.99% (95%CI: 93.11%-99.01%)

Relative Sensitivity

99.59% (95%CI: 98.52%-99.95%)

Relative Specificity

98.92% (95%CI: 97.79%-99.57%)

Accuracy

98.92% (95%CI: 97.67%-99.60%)

Influenza A+B Test:

Description

Test Level

SARS-CoV-2 and Influenza A+B Antigen Combo Rapid Test	RT-PCR	Total
Positive	68	2
Negative	3	485
Total	71	487

95.77% (95%CI: 88.14%-99.12%)

Relative Sensitivity

99.59% (95%CI: 98.52%-99.95%)

Relative Specificity

99.10% (95%CI: 97.92%-99.71%)

Accuracy

98.92% (95%CI: 97.67%-99.60%)

Specificity Testing with Various Viral Strains
The SARS-CoV-2 and Influenza A+B Antigen Combo Rapid Test was tested with the following viral strains. No discernible line at either of the test-line regions was observed at these concentrations listed:

SARS-CoV-2 Test:

Description

Test Level

SARS-CoV-2 and Influenza A+B Antigen Combo Rapid Test	RT-PCR	Total

<tbl_r cells="3" ix="3

